

Dr. DiTolla's

CLINICAL TIPS

PRODUCT *TrollFoil*

CATEGORY *Articulating Paper*

SOURCE *TrollDental*
Newtown Conn.
800-537-8765
trolldentalusa.com

I am unsure when the last significant innovation in articulating paper took place, but I know we have come a long way since the days of typewriter ribbon and cardboard-thick paper. TrollDental hasn't reinvented the articulating wheel with this product, but minor improvements make TrollFoil my favorite articulating paper. First of all, it comes mounted in its own plastic frame, and one less instrument on the bracket table is fine with me. The double-sided foil is only 8 microns thick, and it has no problem marking wet surfaces, dry surfaces or highly polished surfaces, such as cast gold or BruxZir®.

