





See why you should start using TrollFoil today!

Troll®Foil

The small box protects the foil and makes it easy to pick out a single TrollFoil, without getting ink on your fingers.

Individual pieces stacked nicely into a little container that fit nicely in any cabinet, drawer or counter top.



FEATURES & BENEFITS

- Blue or Red foil
- Thickness 8 micron
- Inked on both sides
- Premounted in a frame
- Marks on wet & dry surfaces
- Marks on gold & porcelain
- Timesavings in chair
- Sensitive to pressure



Extremely thin for the highest precision possible

The foil covers half the dental arch. The fact that it is mounted in a semi-rigid frame means that the foil is always flat and stretched.



Distinctive markings on all surfaces

Precise diagnosis, even on wet surfaces, polished surfaces, porcelain and gold. A light tap from the patient and you get the perfect mark.



Watch Webinars & Videos on trolldentalusa.com



REF NUMBER	PRODUCT
12000101	TrollFoil Blue 8 Microns - 100 pcs
12000501	TrollFoil Blue 8 Microns - 500 pcs
12000700	TrollFoil Red 8 Microns - 100 pcs
12000600	TrollFoil Red 8 Microns - 500 pcs